SMART

Kompresörlü Nebulizatör

Model MDS-1000

CE
1984

Türkçe

İçindekiler

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Compressor Nebulizer
Model MDS-1000

Index

Intended Use	18
Who shall use it?	18
Endurance Time	18
How the Nebulizer Works	19
Know Your Unit and Its Accessories	19
Kit Set-Up Diagram	20
Your Device is Ready for First Use	21
Cleaning and Disinfecting	22
Changing the Air Filter	23
Important Safety Instructions	24-25
Troubleshooting	26
Technical Information / Electromagnetic Compatibility (EMC) Statement	27-30
Indications - Contraindications	31
Warranty	32

INTENDED USE

Compressor Nebulizers make the liquid medicine vapor and use a mask. or for inhalation It is the device that will allow it to reach the respiratory tract. Compressor Nebulizer, compressed air and liquid drugs provides available. Compressor Nebulizer means steam generator. Compressor Nebulizers in adults and children worldwide emergency diseases especially of respiratory system diseases It is used. Especially in severe crisis, baby or young child, elderly even your information can be used. Operating the device and putting the mask on the face or close to the respiratory organs is sufficient. Also to take the medicine no action is required. With necessary breathing apparatus assistive devices on sale. Especially the use of drugs to be inhaled used here.

Compressor nebulizers are used to dispense the liquid into particles that the patient can breathe. It is necessary to use compressed air to separate (particles) It is working. Like choose the compressor nebulizer profile it is working.

WHO SHALL USE IT?

- Medical professionals with certificate such as doctors, nurses, therapists.
- At home treatments, patient or caregiver who is under the guidance of a qualified medical specialist.
- In addition the user must be able to understand the content of this

DURABLE PERIOD

Provided the product is used to nebulize 3 times a day for 10 minutes each time at room temperature (23°C), the following durable period is valid. The specified period may vary depending on usage environment.

Durable periods are given below.

Compressor (Main Unit)	10 years
Nebulizer Kit	2 years
Mouthpiece	2 years
Air Tube (PVC, 100 cm)	2 years
Air Filter	60 Applications
Mask (PVC)	2 years

HOW THE NEBULIZER WORKS

This unit is compact, portable and easy to start. It uses a single cylinder piston pump and oil-free lubrication to create Air Pressure.

The air passes through the air filter to the Nebulizer's air compressor unit, to the jet outlet and then the Nebulizer.

,The air flow direction is shown in Figure 1 below.



Figure 1

KNOW YOUR UNIT AND ITS ACCESSORIES

Adult
Mask

Child
Mask


Mouthpiece

Air
Tube


Medication
Cup

Air
Filter

KIT SET-UP DIAGRAM


YOUR DEVICE IS READY FOR FIRST USE

- 1- Make sure the power switch is off.
- 2- Plug in the electrical plug.
3. Remove the inhalation kit cover from the medication cup by gently pulling it up.
4. Put the prescribed amount of medication into the medication cup.
- 5- Make sure that the nebulizer in the medication cup is properly installed.
- 6- Put the medication cup cover back onto the the medication cup.
- 7- Connect the mask or mouthpiece in the nebulizer kit firmly.
- 8- Insert the air tube plug firmly into the air tube connector.
- 9- Turn the power switch to **O**. When the compressor is running, nebulization begins and aerosol is produced. Inhale and exhale while wearing a mask and mouthpiece.

Warning

Do not tilt the nebulizer kit in all directions at an angle greater than 45 degrees. The medication may flow into the mouth or the nebulization may be inefficient.

Correct Use

Incorrect Use

CLEANING AND DISINFECTING

After each use, clean the parts to remove any drug residue. This will prevent ineffective nebulization and reduce the risk of infection.

Washable parts

Nebulizer Kit, Mask (PVC), Mouthpiece

Wash these in hot water or water with vinegar. Rinse thoroughly in clean and hot tap water, tap gently to remove excess water and allow to dry in a clean environment. It is recommended to replace the kit after 100 to 120 applications or approximately 20 washing cycles for a single patient.

Non-Washable Parts

Compressor, Air Tube (PVC)

First of all, make sure that the power plug is not plugged into the outlet. Wipe it with a soft, damp cloth, using water or mild detergent.

Air Filter

Do not wash or clean the air filter. Replace the air filter if it gets wet. Damp air filters can cause clogging.

Disinfecting

Always disinfect the parts after the last treatment of the day is finished. If parts are heavily stained, replace them with new ones.

⚠ Warning : Use of Baffle

Always wash the baffle after each use.

Do not use a brush or needle to clean the baffle.

Make sure you use plenty of water when disinfecting the parts,

Do not wash the baffle with objects other than the appropriate baffle accessories.

Washable Parts

Adult
Mask

Children
Mask

Mouthpiece

Medication
Cup

Baffle


Kit Cover


CHANGING THE AIR FILTER

If the air filter color has changed or has been used for more than 60 applications, replace it with a new one.

1. Remove the air filter cover as shown in the figure.

2. Remove the old air filter with a sharp object such as a toothpick and install the new air filter.

3. Insert the air filter cover into the slot so that it is level with the surface

Notes:

- Only use the air filters that are designed for this device. Do not operate the device without a filter.
- Do not attempt to wash or clean the filter. Replace the air filter if it gets wet. Wet air filters can cause blockages and increase the risk of infection.
- Air filters do not have front or rear side.
- Before installing, verify that the air filters are clean and dust-free.

IMPORTANT SAFETY INSTRUCTIONS

Read the user manual and all other information in the box before using the device.

- It is recommended to use the device for one patient.
- Do not use the device for anything other than nebulization for aerosol treatment. The manufacturer cannot be held responsible for improper use.
- Use only the original accessories specified in the instruction manual for this device.
- Follow the instructions given by your doctor or respiratory therapist regarding the type, dose and regimen of the drug. Check the position of the Nebulization Control ring before use.
- If you feel anything strange during use, stop using the device immediately and consult your doctor.
- Do not use the Nebulizer by placing only water in the medication tank. Keep the device away from the reach of unsupervised children. The device may contain small parts that can be easily swallowed.
- Considering their length, the power cord and the air tubing may cause strangulation hazard.
- Do not store the air tubing while there is moisture or medication remains in it.
- Do not use or store the device where it may be exposed to hazardous fumes or volatile substances.
- Do not use the device where it may be exposed to flammable gases.
- Never use a blanket, towel, etc. to cover the device when in operation.
- Always discard the remaining drug after use and use fresh medicine every time.
- Do not use anesthetic or ventilator breathing circuits.
- Do not block the ventilation slots. Never put the appliance in a place where the ventilation slots may become blocked while operating.
- Keep the parts in a clean place in order to avoid infection.

(Risk of Electric Shock)

- Do not plug in or plug out the device with wet hands.
- The compressor and power cable are not waterproof. Do not spill water or other liquids on these parts. If a liquid is spilled on these parts, immediately unplug the device and wipe off the liquid with a gauze or similar absorbent material.
- Do not immerse the compressor in water or other liquids.
- Do not use or store the appliance in a humid environment.
- Do not use the device with a damaged power cord or plug.
- Keep the power cord away from hot surfaces.

IMPORTANT SAFETY INSTRUCTIONS

(Cleaning and Disinfecting)

Follow the guidelines below when cleaning or disinfecting the parts. Failure to observe these rules may result in damage, inadequate nebulization or infection.

For instructions, see "Cleaning and disinfecting section". Wash and disinfect the nebulizer, mask or mouthpiece before use in the following situations:

- On first use after purchase.
- If the device has not been used for a long period of time.
- If the same device is used by more than one person.

- Make sure that the parts are thoroughly washed and wiped after use and ensure that they are dried and stored in a clean place.
- Be careful not to leave any cleaning solution in the parts.
- Pay extra attention when this device is used by children or people with disabilities, on them or near them.
- Make sure that the parts are connected correctly.
- Ensure that the air filter is properly connected.
- Make sure the air filter is clean. Replace the air filter with a new one, if it has been used for more than 60 applications.
- When using the nebulizer kit, do not tilt it in all angles greater than 45 degrees or shake.
- Do not use or store the device with the air tube bent.
- Only use original nebulization parts, air tubing, air filter and filter cover.
- Do not place more drugs than the medication recommended by your doctor.
- Do not carry or leave the nebulizer kit when the medication cup has medication in it.
- Do not leave the device with children or persons who cannot express their consent.
- This product should not be used in patients who are unconscious or unable to breathe spontaneously.
- Do not expose the device or parts to strong shocks, such as dropping the device.
- Do not distort the baffle or poke the medication cup nozzle with sharp objects such as pins.
- Do not insert fingers or objects into the compressor.
- Do not disassemble or attempt to repair the compressor or power cord.
- Ensure that the baffle is installed correctly before use.


TROUBLESHOOTING

PROBLEM	CAUSE	REMEDY
Nothing happens when the power switch is turned on.	Is the power plug plugged into an electrical outlet?	Control that the plug is inserted in an electrical outlet. Unplug and insert the plug back if necessary.
No nebulization is performed or nebulization rate is low, when the power is on	Is there medication in medication cup?	Add the correct amount of medication to the medication cup.
	Is there too much/little medication in the medication cup?	
	Is the baffle missing and not assembled correctly?	Attach the baffle correctly
	Is the nebulizer kit assembled correctly?	Assemble the nebulizer kit correctly.
	Is the nozzle blocked?	Make sure that the nozzle is free of blockages
	Is the nebulizer kit bent at a sharp angle?	Make sure that the nebulizer kit is not tilted at an angle of more than 45 degrees.
	Is the air tube connected correctly?	Make sure that the air tube is correctly connected to the compressor and nebulizer kit.
	Is the air tube folded or damaged?	Make sure that the air tube is not bent in an angle that it disturbs airflow.
	Is the air tube clogged?	Make sure that the air tube is free of clogging.
	Is the air filter dirty?	Replace the air filter with a new one
The compressor is abnormally loud.	Is the air filter cover attached correctly?	Attach the air filter cover correctly.
The compressor is very hot	Is the compressor covered?	Do not put any covers on it while in operation
	Are the ventilation slots blocked?	Do not block the ventilation slots

TECHNICAL DATA

Voltage / Frequency	AC ~ 230 V / 50-60 HZ	Compressor Pressure Range	2,8 - 3,5 bar
Power / Current	170 VA / 0,7 A	Compressor Flow Rate	8,12 l / dak
Medication Capacity	6 ml	Compressor Outlet Pressure	8 psi
Particle Size	0,4 ile 9 µm	Aerosol Amount	0,08 ml / dak
Mass Median Aerodynamic Diameter (MMAD)	2,79 µm	Operating temperature	0 °C ila +40 °C
Geometrical Standard Deviation	2,57	Humidity	%10 ila %95 relative humidity, concentration
Average Rate of Nebulization	0,07ml	Storage and Transport Temperature	-10 °C ila +70 °C
Noise Level	70 dBA under	Storage and Transport Humidity	%30 by %60 relative humidity, concentration

Result of cascade impactor measurements for particle size

Electromagnetic Compatibility (EMC) Decalaration

This equipment generates, uses and can radiate radio frequency (RF) energy. This equipment may cause electromagnetic interference if it is not installed and used in accordance with the manual.

This device has been tested for Medical Products in accordance with EN-60601-1-2:2015 and has been determined to comply with acceptable limits. These limits indicate that the device provides an acceptable level of protection against electromagnetic interference (EMC) if the device is used as specified in the manual.

This device may be affected by portable and mobile RF communication devices. This equipment must not be stored together with other equipment

See the following tables for more information about this device and EMC.

Guidance and manufacturer's declaration – electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment		
Emission test	Compliance	Electromagnetic environment – guidance
Emissions RF CISPR 11	Group 1	The device is suitable for use in all establishments including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF Emissions CISPR 11	Class B	
Harmonic emissions used for domestic purposes. IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic emissions			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharges IEC 61000-4-2	± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	± 8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ±2 kV Input/output lines: ±1 kV 100 kHz repetition frequency	Power supply lines: ±2 kV Input/output lines: ±1 kV 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surges Line-to-line IEC 61000-4-5	± 0,5 kV, ± 1 kV	± 0,5 kV, ± 1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surges Line-to-ground IEC 61000-4-5	± 0,5 kV, ± 1 kV, ± 2 kV	± 0,5 kV, ± 1 kV, ± 2 kV	
Voltage dips IEC 61000-4-11	0 % U_r ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_r ; 1 cycle and 70% U_r ;25/30 cycle Single phase: 0°	0 % U_T ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_r ;1 cycle and 70% U_r ;25/30 cycle Single phase: 0°	Mains power quality should be that of a typical commercial or hospital environment.
Voltage interruptions IEC 61000-4-11	0% U_r ;250/300 cycle	0% U_r ;250/300 cycle	
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Not: U_r is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic emissions			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	150kHz to 80MHz: 3Vrms 6 Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	150kHz to 80MHz: 3Vrms 6 Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of the device including cables, than the recommended separation distance calculated from the equation applicable to the frequency the transmitter. Recommended separation distance: $d = 0,35\sqrt{P}$ $d = 1,2\sqrt{P}$ 80 MHz to 800 MHz $d = 1,2\sqrt{P}$ 800 MHz to 2.7 GHz $d = 2,3\sqrt{P}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	10V/m 80% Am at 1kHz	10V/m 80% Am at 1kHz	

Note 1: 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- f. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.
- g. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.
- h. Calibration for current injection clamps shall be performed in a 150 Ω system.
- i. If the frequency stepping skips over an ISM or amateur band, as applicable, an additional test frequency shall be used in the ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.
- j. Applicable to ME EQUIPMENT and ME SYSTEMS with RATED input current less than or equal to 16 A / phase and ME EQUIPMENT and ME SYSTEMS with RATED input current greater than 16 A / phase.

Recommended separation distances between portable and mobile RF communications equipment and the device				
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.				
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)			
	150 kHz to 80 MHz (out ISM and amateur radio bands)	150 kHz to 80 MHz (out ISM and amateur radio bands)	80 MHz to 800 MHz	800 MHz to 2.7 GHz
	$d = 1,2\sqrt{P}$	$d = 0,6\sqrt{P}$	$d = 1,2\sqrt{P}$	$d = 2,3\sqrt{P}$
0.01	0.12	0.06	0.12	0.23
0.1	0.38	0.19	0.38	0.73
1	1.2	0.6	1.2	2.3
10	3.8	1.9	3.8	7.3
100	12	6	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.				
Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.				
Note 2: The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.				
Note 3: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

Guidance and manufacturer's declaration – electromagnetic emissions							
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.							
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communica- tions equipments	Test Frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modulation b) (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
	385	380-390	TETRA 400	Pulse modulation b) 18 Hz	1.8	0.3	27
	450	380-390	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	2	0.3	28
	710	704-787	LTE Band 13.17	Pulse modulation b) 217 Hz	0.2	0.3	9
	745						
	780						
	810	800-960	GSM 800/900, TETRA 800, iDEN 820, LTE Band 5	Pulse modulation b) 18 Hz	2	0.3	28
	870						
	930						
	1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation b) 217 Hz	2	0.3	28
	1845						
	1970						
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0.3	28
	5240	5100-5800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	0.2	0.3	9
	5240						
	5785						

Note: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.							
d. For some services, only the uplink frequencies are included.							
e. The carrier shall be modulated using a 50 % duty cycle square wave signal.							
f. As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.							
The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation: $E = 6/d\sqrt{P}$ Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.							

Classification According to EN 60601-1:2006/A1:2013

☐ Category II (Double Insulation)

 Type B Equipment

IP21 Degree of Protection Against Dripping Water

Indications

Indications

This type of drug intake system is recommended for drugs to be inhaled directly as aerosol, but it is suitable for drugs that can be converted into aerosol form.

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Compressor Nebulizers can be used to treat the following conditions:

- Asthma
- Chronic Obstructive Pulmonary Diseases
- Cystic Fibrosis
- Respiratory Infections
- Room Humidification
- Humidification of breathing air Drug Nebulization
- Specific and Non-specific lung diseases
- Bronchial Asthma
- Acute and Chronic Bronchitis
- Pneumonia

Contraindications

Sensitivity and cardiotoxicity to propellants are considered as a contraindication to the use of this device.

Inhaler should not be used for patients who do not provide adequate inhalation during breathing. It should not be used in newborns and premature children who are unconscious

The device has no direct side effects. The important thing is to use medication according to the prescription given by the doctor.

 CE Mark, 93/42/EEC Directive suitable product(2292)	 Store in a dry place.
 B Type Equipment	SN Serial number
 Category II Medical Device (Double Insulation)	IP21 Degree of Protection Against Dripping Water
 See user manual.	 Manufacturer Company Information
 Attention ! Some for this device There are warnings and precautions.	AC  Alternative current
 2011 / 65 / EU ve 2012 / 19 / EU According to the directive, this device is and recycling of electronic equipment should be taken to the points.	

WARRANTY

Thank you for buying the product, MEDİSEC SMART MDS-1000. This product is manufactured with great care and superior quality materials. It is designed to meet your requirements as long as it is used correctly and maintained in accordance with the instructions in the manual.

MEDİSEC provides a 2-year warranty for this product, which starts on the date of purchase. The production, workmanship and material quality of this product is guaranteed by MEDİSEC. During this warranty, MEDİSEC will not demand any costs for parts, labor, repair or replacement of problematic products or parts.

Conditions not covered by the warranty:

- Shipping risks and costs.
- Repair costs incurred by unauthorized persons and / or defects resulting from repair by unauthorized persons.
- Periodic maintenance and controls.
- Failure or wear of the optional parts or parts other than the main unit, apart from the ones within scope of the warranty which are expressly stated above.

If you need a service under warranty, contact the dealer from whom you purchased the product or an authorized MEDİSEC dealer. For the address, refer to the product box, booklets, or the retailer.

Repairs or replacements under the warranty do not result in an extension or renewal of the warranty period. Warranty service is provided if the product is shipped with the original invoice/receipt given to the consumer from the seller. MEDİSEC reserves the right to refuse warranty service if the information provided is not clear.

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